


Scan Clients Anywhere with
ZYTO REMOTE

ZYTO Remote allows you to provide the same valuable information and services you offer without being limited by your physical location.

- Assist clients that are not able to come to you
- Maintain relationships with clients that have moved away
- Run scans on behalf of other wellness professionals
- Expand your business without having to expand office space



Available in the Balance, Select, Elite, & EVOX

HOW IT WORKS



1. Purchase a Remote Hand Cradle.

2. We can ship it directly to your client, or you can receive and ship it yourself.



3. Instruct your client to download and install the Remote app at zyto.com/remote.

4. Click the Remote icon located in the upper-right corner of your software. Then click Connect.




5. Share the remote code with your client.
(Make sure their Hand Cradle is plugged in and their Remote app is opened.)



6. After the client enters the code, instruct them to click Connect and place their hand on the Hand Cradle.

7. Run a scan for the client as you normally would.



Remote session recommendations

- High-speed Internet connection
- Windows 10 operating system
(Note: ZYTO Remote is not Mac-compatible.)
- Videoconferencing/screensharing software (Zoom, Join.Me, Skype, etc.)
- Phone or instant messaging software (WhatsApp, Slack, Skype, etc.)

For current pricing & to place a Remote Hand Cradle order, email sales@zyto.com or call us at (801) 224-7199, option 2.